

Since 1954, SPEX CertiPrep is the industry leader in the CRM marketplace meeting the needs of laboratories worldwide with innovation and research. Accredited by A2LA to ISO/IEC 17025:2017 & ISO 17034:2016. Certified by DQS to ISO 9001:2015.

Ask Us About Our 10 ppm Singles

Claritas PPT® Grade Certified Reference Materials are designed for ICP-MS and can also be used in ICP analysis. We are now offering single-element solutions at 10 ppm with accuracy within 2% of uncertainty. They are made using ultra-high purity acids, the highest grade starting materials and a high purity water to minimize contaminants. For additional product information, please contact us at [732.549.7144](tel:732.549.7144) or via email at CRMSales@spex.com.

Element	Concentration	Volume	Matrix	Part #
Aluminum	10 µg/mL	30 mL	2% HNO ₃	CLAL2-1AM
Aluminum	10 µg/mL	125 mL	2% HNO ₃	CLAL2-1AY
Arsenic	10 µg/mL	30 mL	2% HNO ₃	CLAS2-1AM
Arsenic	10 µg/mL	125 mL	2% HNO ₃	CLAS2-1AY
Cadmium	10 µg/mL	30 mL	2% HNO ₃	CLCD2-1AM
Cadmium	10 µg/mL	125 mL	2% HNO ₃	CLCD2-1AY
Chromium	10 µg/mL	30 mL	2% HNO ₃	CLCR2-1AM
Chromium	10 µg/mL	125 mL	2% HNO ₃	CLCR2-1AY
Copper	10 µg/mL	30 mL	2% HNO ₃	CLCU2-1AM
Copper	10 µg/mL	125 mL	2% HNO ₃	CLCU2-1AY
Gold	10 µg/mL	30 mL	2% HCl	CLAU1-1AM
Gold	10 µg/mL	125 mL	2% HCl	CLAU1-1AY
Iron	10 µg/mL	30 mL	2% HNO ₃	CLFE2-1AM
Iron	10 µg/mL	125 mL	2% HNO ₃	CLFE2-1AY
Lead	10 µg/mL	30 mL	2% HNO ₃	CLPB2-1AM
Lead	10 µg/mL	125 mL	2% HNO ₃	CLPB2-1AY
Manganese	10 µg/mL	30 mL	2% HNO ₃	CLMN2-1AM
Manganese	10 µg/mL	125 mL	2% HNO ₃	CLMN2-1AY

10 ppm Single-Element Standards

Element	Concentration	Volume	Matrix	Part #
Nickel	10 µg/mL	30 mL	2% HNO ₃	CLNI2-1AM
Nickel	10 µg/mL	125 mL	2% HNO ₃	CLNI2-1AY
Palladium	10 µg/mL	30 mL	2% HCl	CLPD1-1AM
Palladium	10 µg/mL	125 mL	2% HCl	CLPD1-1AY
Platinum	10 µg/mL	30 mL	2% HCl	CLPT1-1AM
Platinum	10 µg/mL	125 mL	2% HCl	CLPT1-1AY
Zinc	10 µg/mL	30 mL	2% HNO ₃	CLZN2-1AM
Zinc	10 µg/mL	125 mL	2% HNO ₃	CLZN2-1AY

CAN'T FIND THE STANDARDS YOU ARE LOOKING FOR?

SPEX CertiPrep can make custom standards to meet your exact needs. Contact us for more information.

Contact Us

Phone: 800.LAB.SPEX • 732.549.7144 • Fax: 732.603.9647
CRMSales@spex.com • spexcertiprep.com

© 2020 SPEX CertiPrep. All Rights Reserved.

CONNECT WITH US

