

Since 1954, SPEX CertiPrep is the industry leader in the CRM marketplace meeting the needs of laboratories worldwide with innovation and research. Accredited by A2LA to ISO/IEC 17025:2017 & ISO 17034:2016. Certified by DQS to ISO 9001:2015.

Inorganic Product Part Number Formulation

Interested in how SPEX CertiPrep develops our part number for our Certified Reference Materials? Information below will spell out the method to our madness!



Product Type: **PLAL2-2X**

AS: Anion Standards
CL: Claritas
CS: Cation Standards

IC: Ion Chromatography
MM: Matrix Modifier
ORG: Organometallic Oils

PL: Plasma
QC: Quality Control



Element: **PLAL2-2X**

- References elements from the periodic table
- Aluminum to Zinc
- Anions and Cations



Matrix: **PLAL2-2X**

1: 2-5% HCl	4: 10% HNO3	7: HNO2/tr. HNO3/0.6%	X Xylene
2: 2-5% HNO3	5: 20% HCl (packaged in PTFE)	8: Base Oil	
3: 10% HCl	6: 20% HNO3	9: H2O	



Concentration: **PLAL2-2X**

1A: 10 ppm	1: 100 ppm
2: 1,000 ppm	3: 10,000 ppm



Concentration: **PLAL2-2X**

Z: 50 mL	Y: 100 mL
T: 250 mL	X: 500 mL

Contact Us

Phone: 800.LAB.SPEX • 732.549.7144 • Fax: 732.603.9647
CRMSales@spex.com • spexcertiprep.com

© 2020 SPEX CertiPrep. All Rights Reserved.

CONNECT WITH US

